

ENPEP / Aminopeptidase A Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C12167 - LSBio	
CatalogID:	LS-C12167
Target:	glutamyl aminopeptidase (aminopeptidase A) (ENPEP)
Synonyms:	ENPEP Antibody, APA Antibody, Aminopeptidase A Antibody, AP-A Antibody, CD249 antigen Antibody, Differentiation antigen gp160 Antibody, EAP Antibody, Glutamyl aminopeptidase Antibody, gp160 Antibody, CD249 Antibody
Family / Subfamily:	Exopeptidase / Metallopeptidase M1
Host	ENPEP antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ENPEP / Aminopeptidase A antibody was raised against Mouse
Specificity:	Recognizes BP-1/Ly-51/6C3 antigen, (aminopeptidase A, E. C.3.4.11.7, glutamylaminopeptidase).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/24/2014 © 2014 LifeSpan BioSciences	