

MYC / c-Myc Mouse Monoclonal Antibody - LS-C122773 - LSBio

CatalogID:	LS-C122773
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide corresponding to myc sequence.
Specificity:	Recognizes Myc tag.
Purification:	Protein A purified
Presentation:	10 mM phosphate, pH 7.4, 150 mM sodium chloride, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C.
Usage Summary:	Suitable for use in Western Blot and Immunohistochemistry (frozen).
Uses:	IHC - Paraffin, IHC - Frozen, Western blot (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences