

CD69 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C123096 - LSBio

CatalogID:	LS-C123096
Target:	CD69 molecule
Synonyms:	CD69 Antibody, Activation inducer molecule Antibody, AIM Antibody, CD69 antigen Antibody, CLEC2C Antibody, EA1 Antibody, Early activation antigen CD69 Antibody, BL-AC/P26 Antibody, CD69 molecule Antibody, gp32/28 Antibody, MLR-3 Antibody
Host	CD69 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD69 antibody was raised against Human
Immunogen:	CD69 antibody was raised against activated human B-cells.
Specificity:	Recognizes human CD69. Species cross-reactivity: baboon, chimpanzee, cynomolgus monkey, rhesus monkey, macaque.
Reactivity:	Human, Chimpanzee, Baboon, Monkey, Cynomolgus, Macaque
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.2% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. FITC conjugates are sensitive to light.
Usage Summary:	Suitable for use in Flow Cytometry, Immunofluorescence, and Immunohistochemistry (frozen). Flow Cytometry: Neat; 20 ul labels 10 ⁶ cells in 100ul.
Uses:	IHC - Frozen, Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences