

THY1 / CD90 Mouse anti-Rat Monoclonal Antibody - LS-C123118 - LSBio

CatalogID:	LS-C123118
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	THY1 / CD90 antibody was raised against Rat
Antigen Type:	Purified protein
Immunogen:	THY1 / CD90 antibody was raised against rat Thy1 antigen.
Specificity:	Recognizes rat CD90. Species cross-reactivity: Rabbit, Mouse, Guinea Pig.
Reactivity:	Rat, Mouse, Guinea pig, Rabbit
Purification:	Protein G purified
Presentation:	PBS, pH 7.2. No preservative added
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C.
Usage Summary:	Suitable for use in Flow Cytometry, Immunohistochemistry (frozen), Immunoprecipitation, and Western Blot. Flow Cytometry: 1:40-1:80; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry (1:40 - 1:80) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences