

CBLL1 / HAKAI Rabbit anti-Human Polyclonal Antibody - LS-C123676 - LSBio

CatalogID:	LS-C123676
Target:	Cbl proto-oncogene-like 1, E3 ubiquitin protein ligase (CBLL1)
Synonyms:	CBLL1 Antibody, C-Cbl-like protein 1 Antibody, E-cadherin binding protein E7 Antibody, HAKAI Antibody, RING finger protein 188 Antibody, RNF188 Antibody
Host	CBLL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	CBLL1 / HAKAI antibody was raised against Human
Immunogen:	CBLL1 / HAKAI antibody was raised against human Hakai.
Specificity:	Recognizes human Hakai.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 0.1% BSA, 0.02% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months.
Usage Summary:	Suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot. Immunohistochemistry: 1:50. Immunofluorescence: 1:50. Western Blot: 1:50.
Uses:	IHC (1:50), Immunofluorescence (1:50), Western blot (1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences