

MYC / c-Myc Rabbit anti-Human Polyclonal Antibody - LS-C124194 - LSBio

CatalogID:	LS-C124194
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic Non-phosphopeptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic non-phosphopeptide corresponding to human Myc around the phosphorylation site of threonine 358 (R-R-TP-H-N).
Specificity:	Recognizes human c-Myc. Species cross-reactivity: mouse, rat.
Reactivity:	Human, Chimpanzee, Gorilla, Orangutan, Monkey, Mouse, Rat, Bovine, Dog, Hamster, Pig, Rabbit, Chicken, Xenopus
Purification:	Immunoaffinity purified
Presentation:	PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months.
Usage Summary:	Suitable for use in Western Blot and Immunohistochemistry. Western Blot: 1:500-1:1000. Immunohistochemistry: 1:50-1:100.
Uses:	IHC (1:50 - 1:100), Western blot (1:500 - 1:1000) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences