

TLR2 Mouse anti-Mouse Monoclonal (FITC) Antibody - LS-C124666 - LSBio	
CatalogID:	LS-C124666
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR2 antibody was raised against Mouse
Immunogen:	TLR2 antibody was raised against mouse TLR2.
Specificity:	Recognizes mouse TLR2. Species cross-reactivity: human.
Reactivity:	Mouse, Human
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. FITC conjugates are sensitive to light.
Usage Summary:	Suitable for use in Immunohistochemistry and Flow Cytometry. Immunohistochemistry (frozen): 1:50. Flow Cytometry: 1:50.
Uses:	IHC - Paraffin, IHC - Frozen (1:50), Flow Cytometry (1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	