

RP105 / CD180 Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C125327 - LSBio

CatalogID:	LS-C125327
Target:	CD180 molecule
Synonyms:	CD180 Antibody, CD180 antigen Antibody, Lymphocyte antigen 64 Antibody, Ly78 Antibody, LY64 Antibody, CD180 molecule Antibody, RP105 Antibody
Host	CD180 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RP105 / CD180 antibody was raised against Mouse
Immunogen:	RP105 / CD180 antibody was raised against mouse CD 180.
Specificity:	Recognizes mouse RP 105.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 1:50.
Uses:	Flow Cytometry (1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences