

TLR1 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C126372 - LSBio	
CatalogID:	LS-C126372
Target:	toll-like receptor 1 (TLR1)
Synonyms:	TLR1 Antibody, CD281 antigen Antibody, KIAA0012 Antibody, Rsc786 Antibody, Toll-like receptor 1 Antibody, TIL. LPRS5 Antibody, TIL Antibody, CD281 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR1 antibody was raised against Human
Immunogen:	TLR1 antibody was raised against human TLR1.
Specificity:	Recognizes human TLR1. Does not recognizes TLR2, TLR4 (CD284) or TLR4 in association with MD2.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 1:50.
Uses:	Flow Cytometry (1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	