

TLR6 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C126387 - LSBio

CatalogID:	LS-C126387
Target:	toll-like receptor 6 (TLR6)
Synonyms:	TLR6 Antibody, CD286 antigen Antibody, CD286 Antibody, Toll-like receptor 6 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR6 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR6 antibody was raised against Human
Immunogen:	TLR6 antibody was raised against human TLR6.
Specificity:	Recognizes human TLR6.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Immunofluorescence, and Flow Cytometry. Immunofluorescence: 1:50. Flow Cytometry: 1:50.
Uses:	Immunofluorescence (1:50), Flow Cytometry (1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences