

**V5 Tag Mouse Monoclonal (FITC) Antibody - LS-C126452 - LSBio**

<b>CatalogID:</b>	LS-C126452
<b>Target:</b>	V5 Tag
<b>Host</b>	V5 Tag antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	V5 Tag antibody was raised against recombinant peptide corresponding to the PK tag fused to a class II restricted T helper cell epitope.
<b>Specificity:</b>	Recognizes a small epitope, termed Pk, present on the P/V proteins of the paramyxovirus, SV5.
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 1% BSA
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
<b>Usage Summary:</b>	Suitable for use in Immunofluorescence and Suitable for use in Flow cytometry: Neat-1:5.
<b>Uses:</b>	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences