

CCR5 Rat anti-Human Monoclonal (PE) Antibody - LS-C127289 - LSBio	
<b>CatalogID:</b>	LS-C127289
<b>Target:</b>	chemokine (C-C motif) receptor 5 (gene/pseudogene) (CCR5)
<b>Synonyms:</b>	CCR5 Antibody, C-C chemokine receptor 5 Antibody, CD195 antigen Antibody, CD195 Antibody, C-C CKR-5 Antibody, Chemokine c-c motif receptor 5 Antibody, Chemr13 Antibody, CMKBR5 Antibody, CKR5 Antibody, Cc chemokine receptor 5 Antibody, CCCKR5 Antibody, CCR-5 Antibody, Ccr5 receptor Antibody, Chemokine receptor CCR5 Antibody, HIV-1 fusion coreceptor Antibody, IDDM22 Antibody, C-C chemokine receptor type 5 Antibody, CC-CKR-5 Antibody, CKR-5 Antibody
<b>Family / Subfamily:</b>	GPCR / Chemokine
<b>Host</b>	CCR5 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CCR5 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CCR5 antibody was raised against CHO cells transfected with human CCR5.
<b>Specificity:</b>	Recognizes human CCR5 at ~45kD.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Reconstitution:</b>	Sterile distilled water or PBS
<b>Presentation:</b>	Lyophilized from PBS, 1% BSA, 0.09% sodium azide
<b>Recommended Storage:</b>	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
<b>Usage Summary:</b>	Suitable for use in Flow cytometry: Neat-1:5; 10 ul labels 10 <sup>6</sup> cells in 100ul.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/24/2014 © 2014 LifeSpan BioSciences	