

CD4 Rat anti-Dog Monoclonal (FITC) Antibody - LS-C127352 - LSBio

CatalogID:	LS-C127352
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD4 antibody was raised against Dog
Antigen Type:	Cells
Immunogen:	CD4 antibody was raised against canine Concanavalin A activated T-cell blasts.
Specificity:	Recognizes canine CD4.
Reactivity:	Dog
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences