

CD4a Mouse anti-Pig Monoclonal (FITC) Antibody - LS-C127362 - LSBio

CatalogID:	LS-C127362
Target:	CD4a
Host	CD4a antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD4a antibody was raised against Pig
Immunogen:	CD4a antibody was raised against leukocytes isolated from porcine gut lamina propria.
Specificity:	Recognizes porcine CD4a.
Reactivity:	Pig
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences