

CD4 Mouse anti-Pig Monoclonal (PE) Antibody - LS-C127363 - LSBio

CatalogID:	LS-C127363
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD4 antibody was raised against Pig
Immunogen:	CD4 antibody was raised against leukocytes isolated from porcine gut lamina propria.
Specificity:	Recognizes porcine CD4a.
Reactivity:	Pig
Purification:	Protein G purified
Reconstitution:	Sterile buffer or distilled water
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences