

CD5 Rat anti-Dog Monoclonal (PE) Antibody - LS-C127375 - LSBio

CatalogID:	LS-C127375
Target:	CD5 molecule
Synonyms:	CD5 Antibody, CD5 antigen Antibody, CD5 antigen (p56-62) Antibody, CD5 molecule Antibody, Lymphocyte antigen T1/Leu-1 Antibody, LEU1 Antibody, T1 Antibody
Host	CD5 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD5 antibody was raised against Dog
Immunogen:	CD5 antibody was raised against concanavalin A activated canine peripheral blood cells.
Specificity:	Recognizes canine CD5, a 67kD cell surface glycoprotein expressed on the surface of T-cells.
Reactivity:	Dog
Purification:	Protein G purified
Reconstitution:	Sterile buffer or distilled water
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences