

**ITGAM / CD11b Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C127440 - LSBio**

<b>CatalogID:</b>	LS-C127440
<b>Target:</b>	integrin, alpha M (complement component 3 receptor 3 subunit) (ITGAM)
<b>Synonyms:</b>	ITGAM Antibody, Antigen CD11b (p170) Antibody, CD11B Antibody, CD11b antigen Antibody, CR3A Antibody, Integrin alpha-M Antibody, MAC-1 Antibody, MO1A Antibody, Neutrophil adherence receptor Antibody, MAC1A Antibody, SLEB6 Antibody, CR-3 alpha chain Antibody
<b>Family / Subfamily:</b>	Integrin / not assigned-Integrin
<b>Host</b>	ITGAM antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ITGAM / CD11b antibody was raised against Mouse
<b>Immunogen:</b>	ITGAM / CD11b antibody was raised against thioglycollate-elicited Peritoneal Macrophages (TPM).
<b>Specificity:</b>	Recognizes mouse CD11b. Species cross-reactivity: human.
<b>Reactivity:</b>	Mouse, Human
<b>Purification:</b>	Ion exchange chromatography
<b>Presentation:</b>	PBS, pH 7.2, 1% BSA, 0.09% sodium azide
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
<b>Usage Summary:</b>	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 <sup>6</sup> cells in 100ul.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences