

CD22 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C127539 - LSBio	
CatalogID:	LS-C127539
Target:	CD22 molecule
Synonyms:	CD22 Antibody, CD22 molecule Antibody, CD22 antigen Antibody, SIGLEC2 Antibody, SIGLEC-2 Antibody, T-cell surface antigen Leu-14 Antibody, B-cell receptor CD22 Antibody, BL-CAM Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD22 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD22 antibody was raised against Mouse
Immunogen:	CD22 antibody was raised against murine CD22.
Specificity:	Recognizes domain 2 of the murine cell surface glycoprotein CD22.
Reactivity:	Mouse
Purification:	Protein G purified
Reconstitution:	Sterile buffer or distilled water
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat-1:5; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	