

CD34 Mouse anti-Dog Monoclonal (PE) Antibody - LS-C127656 - LSBio

CatalogID:	LS-C127656
Target:	CD34 molecule
Synonyms:	CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody
Host	CD34 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD34 antibody was raised against Dog
Immunogen:	CD34 antibody was raised against canine CD34 fusion protein.
Specificity:	Recognizes canine CD34, a cell surface antigen of approximately 110kD expressed by endothelial cells and hematopoietic stem cells.
Reactivity:	Dog
Purification:	Protein G purified
Reconstitution:	Sterile buffer or distilled water
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences