

CD55 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C127776 - LSBio	
CatalogID:	LS-C127776
Target:	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
Synonyms:	CD55 Antibody, CD55 antigen Antibody, Decay-accelerating factor Antibody, TC Antibody, CR Antibody, CROM Antibody, DAF Antibody
Host	CD55 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD55 antibody was raised against Human
Immunogen:	CD55 antibody was raised against human CD55.
Specificity:	Recognizes human CD55.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	