

SELL / L-Selectin / CD62L Mouse anti-Bovine Monoclonal (FITC) Antibody - LS-C127804 - LSBio	
CatalogID:	LS-C127804
Target:	selectin L (SELL)
Synonyms:	SELL Antibody, CD62L Antibody, CD62L antigen Antibody, LAM-1 Antibody, LAM1 Antibody, LECAM1 Antibody, LEU8 Antibody, Leukocyte adhesion molecule 1 Antibody, Lyam-1 Antibody, LSEL Antibody, gp90-MEL Antibody, HLHRc Antibody, PLNHR Antibody, L-selectin Antibody, Leu-8 Antibody, Selectin L Antibody, LNHR Antibody, LYAM1 Antibody, Lymph node homing receptor Antibody, Lymphocyte adhesion molecule 1 Antibody, TQ1 Antibody, PIn homing receptor Antibody
Host	SELL antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SELL / L-Selectin / CD62L antibody was raised against Bovine
Immunogen:	SELL / L-Selectin / CD62L antibody was raised against bovine CD62L.
Specificity:	Recognizes bovine CD62L.
Reactivity:	Bovine
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	