

**CD71 / Transferrin Receptor Mouse anti-Human Monoclonal (PE) (MEM75) Antibody - LS-C127833
- LSBio**

CatalogID:	LS-C127833
Target:	transferrin receptor (TFRC)
Synonyms:	TFRC Antibody, CD71 Antibody, TFR1 Antibody, TRFR Antibody, TFR Antibody, TR Antibody, CD71 antigen Antibody, p90 Antibody, T9 Antibody, Transferrin receptor Antibody, Transferrin receptor protein 1 Antibody
Family / Subfamily:	Protease / Metallopeptidase M28B
Host	TFRC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM75
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD71 / Transferrin Receptor antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD71 / Transferrin Receptor antibody was raised against pre-B cell line NALM-6.
Specificity:	Recognizes human CD71. Does not block binding of transferrin to the receptor.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.4, BSA, 15 mM sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 20 ul labels human blood cells in 100 ul whole blood.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences