

CD163 Mouse anti-Rat Monoclonal (FITC) Antibody - LS-C128001 - LSBio

CatalogID:	LS-C128001
Target:	CD163 molecule
Synonyms:	CD163 Antibody, CD163 antigen Antibody, Macrophage-associated antigen Antibody, M130 Antibody, CD163 molecule Antibody, Hemoglobin scavenger receptor Antibody, MM130 Antibody
Host	CD163 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD163 antibody was raised against Rat
Immunogen:	CD163 antibody was raised against rat spleen cell homogenate.
Specificity:	Recognizes rat CD163. Does not react with human.
Reactivity:	Rat
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 1% BSA
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences