

ITGAX / CD11c Mouse anti-Human Monoclonal (PE) Antibody - LS-C12949 - LSBio

CatalogID:	LS-C12949
Target:	integrin, alpha X (complement component 3 receptor 4 subunit) (ITGAX)
Synonyms:	ITGAX Antibody, CD11c antigen Antibody, Integrin alpha-X Antibody, Leu M5, alpha subunit Antibody, p150 95 integrin alpha chain Antibody, SLEB6 Antibody, CD11C Antibody, Leu M5 Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGAX antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ITGAX / CD11c antibody was raised against Human
Specificity:	Recognizes human CD11c.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Flow Cytometry: 20 ul/1E6 cells.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences