

Human IgD Goat anti-Human Polyclonal (PE) Antibody - LS-C129915 - LSBio	
CatalogID:	LS-C129915
Target:	Human IgD
Host	Human IgD antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Human IgD antibody was raised against Human
Immunogen:	Human IgD antibody was raised against human IgD.
Specificity:	Recognizes human IgD.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in ELISA and Flow Cytometry.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	250 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	