

Hemagglutinin / HA Tag Chicken Polyclonal (FITC) Antibody - LS-C130008 - LSBio

CatalogID:	LS-C130008
Target:	HA Tag / Hemagglutinin Tag
Host	Chicken
Clonality:	Polyclonal
Isotype:	IgY
Conjugations:	Fluorescein (FITC)
Immunogen:	Influenza hemagglutinin-HA-epitope (YPYDVPDYA) (KLH).
Specificity:	Recognizes Hemagglutinin Epitope Tag.
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% BSA, 0.05% sodium azide, 50% glycerol
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Usage Summary:	Suitable for use in ELISA and Western Blot.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences