

IL4 Mouse anti-Bovine Monoclonal (PE) Antibody - LS-C130198 - LSBio	
CatalogID:	LS-C130198
Target:	interleukin 4 (IL4)
Synonyms:	IL4 Antibody, B cell growth factor 1 Antibody, B-cell stimulatory factor 1 Antibody, B_cell stimulatory factor 1 Antibody, BCGF1 Antibody, BSF-1 Antibody, Binetrakin Antibody, BCGF-1 Antibody, Interleukin 4 Antibody, Interleukin-4 Antibody, Pitrakinra Antibody, BSF1 Antibody, IL-4 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL4 antibody was raised against Bovine
Immunogen:	IL4 antibody was raised against bovine Interleukin 4.
Specificity:	Recognizes bovine interleukin 4. Species cross-reactivity: Dog, Pig, Sheep, Mustelid, Goat.
Reactivity:	Bovine, Dog, Goat, Pig, Sheep, Mustelid
Purification:	Protein G purified
Reconstitution:	Restore with 1ml distilled water.
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	