

KLH Chicken anti-Keyhole Limpet Polyclonal (FITC) Antibody - LS-C130346 - LSBio	
CatalogID:	LS-C130346
Target:	KLH / Keyhole Limpet Hemocyanin
Host	KLH / Keyhole Limpet Hemocyanin antibody was produced in Chicken
Clonality:	Polyclonal
Isotype:	IgY
Conjugations:	Fluorescein (FITC)
Immunogen Species:	KLH antibody was raised against Keyhole Limpet
Immunogen:	KLH antibody was raised against purified Keyhole Limpet Hemocyanin (KLH).
Specificity:	Recognizes Keyhole Limpet Hemocyanin.
Reactivity:	Keyhole Limpet
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.01% Thimerosal, 50% glycerol
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Usage Summary:	Suitable for use in ELISA.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	