

MSR1 / CD204 Rat anti-Mouse Monoclonal (Extracellular Domain) (PE) (268318) Antibody - LS-C130609 - LSBio	
CatalogID:	LS-C130609
Target:	macrophage scavenger receptor 1 (MSR1)
Synonyms:	MSR1 Antibody, CD204 antigen Antibody, SCARA1 Antibody, PhSR1 Antibody, SR-A Antibody, CD204 Antibody, PhSR2 Antibody, SRA Antibody
Host	MSR1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	268318
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MSR1 / CD204 antibody was raised against Mouse
Immunogen:	MSR1 / CD204 antibody was raised against nS0-derived, recombinant mouse Scavenger Receptor class A isoform I (rmSR-AI; aa83-458) extracellular domain.
Specificity:	Recognizes mouse MSR1.
Epitope:	Extracellular Domain
Reactivity:	Mouse
Purification:	Immunoaffinity purified
Presentation:	Saline, 0.5% BSA, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	