

CD20 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C13088 - LSBio	
CatalogID:	LS-C13088
Target:	membrane-spanning 4-domains, subfamily A, member 1 (MS4A1)
Synonyms:	MS4A1 Antibody, B-lymphocyte antigen CD20 Antibody, B1 Antibody, Bp35 Antibody, CD20 Antibody, CD20 antigen Antibody, CD20 receptor Antibody, LEU-16 Antibody, CVID5 Antibody
Host	MS4A1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD20 antibody was raised against Human
Immunogen:	CD20 antibody was raised against human CD20.
Specificity:	Recognizes human CD20.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry and immunohistochemistry. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	