

NKG2A/C/E Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C131065 - LSBio

CatalogID:	LS-C131065
Target:	NKG2A/C/E
Host	NKG2A/C/E antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NKG2A/C/E antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	NKG2A/C/E antibody was raised against CHO transfected cells expressing the B6 allele of NKG2A.
Specificity:	Recognizes mouse NKG2A, NKG2C and NKG2E, which are isoforms of the NKG2 lectin-like family. Reported to block ligand binding to the receptor.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat-1:5; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences