

**Protein L Chicken anti-Peptostreptococcus magnus Polyclonal (FITC) Antibody - LS-C131247 - LSBio**

<b>CatalogID:</b>	LS-C131247
<b>Target:</b>	Protein L
<b>Host</b>	Protein L antibody was produced in Chicken
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgY
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Protein L antibody was raised against Peptostreptococcus magnus
<b>Immunogen:</b>	Protein L antibody was raised against purified recombinant protein L (cell surface protein from Peptostreptococcus magnus).
<b>Specificity:</b>	Recognizes Protein L.
<b>Reactivity:</b>	Peptostreptococcus magnus
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.05% sodium azide, 50% glycerol
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Usage Summary:</b>	Suitable for use in ELISA.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences