

FCER2 / CD23 Mouse anti-Human Monoclonal (PE) Antibody - LS-C13133 - LSBio

CatalogID:	LS-C13133
Target:	Fc fragment of IgE, low affinity II, receptor for (CD23) (FCER2)
Synonyms:	FCER2 Antibody, CD23 antigen Antibody, CD23A Antibody, CLEC4J Antibody, CD23 Antibody, Fc-epsilon-RII Antibody, IGEBF Antibody, Lymphocyte IgE receptor Antibody, BLAST-2 Antibody, FCE2 Antibody
Host	FCER2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG3
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FCER2 / CD23 antibody was raised against Human
Specificity:	Recognizes human CD23.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 5 ul/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences