

Phosphotyrosine Mouse Monoclonal (PE) Antibody - LS-C131435 - LSBio

CatalogID:	LS-C131435
Target:	Phosphotyrosine
Host	Phosphotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Phycoerythrin (PE)
Antigen Type:	Hapten
Immunogen:	Phosphotyrosine antibody was raised against phosphotyrosine.
Specificity:	Recognizes Phosphotyrosine.
Purification:	Purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences