

SIGLEC1 / CD169 / Sialoadhesin Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C132105 - LSBio

CatalogID:	LS-C132105
Target:	sialic acid binding Ig-like lectin 1, sialoadhesin (SIGLEC1)
Synonyms:	SIGLEC1 Antibody, CD169 Antibody, CD169 antigen Antibody, SIGLEC-1 Antibody, Sialoadhesin Antibody, SN Antibody, DJ1009E24.1 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	SIGLEC1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SIGLEC1 / CD169 / Sialoadhesin antibody was raised against Mouse
Immunogen:	SIGLEC1 / CD169 / Sialoadhesin antibody was raised against stromal (reticular) elements from spleen.
Specificity:	Recognizes mouse Siglec-1. Does not react with human or rat.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.09% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences