

CD30 Hamster anti-Mouse Monoclonal (Extracellular Domain) (FITC) Antibody - LS-C13243 - LSBio	
CatalogID:	LS-C13243
Target:	tumor necrosis factor receptor superfamily, member 8 (TNFRSF8)
Synonyms:	TNFRSF8 Antibody, CD30L receptor Antibody, CD30 Antibody, Cytokine receptor CD30 Antibody, D1S166E Antibody, Ki-1 Antibody, Ki-1 antigen Antibody, CD30 antigen Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF8 antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD30 antibody was raised against Mouse
Immunogen:	CD30 antibody was raised against recombinant mouse CD30 extracellular domain-mouse IgG1 fusion protein.
Specificity:	Recognizes murine CD30, Mr 52kD (core peptide).
Epitope:	Extracellular Domain
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. In vitro functional studies. Flow Cytometry: 3 ug/1E6 cells.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	