

WC1 Surface Antigen Mouse anti-Bovine Monoclonal (FITC) Antibody - LS-C132931 - LSBio

CatalogID:	LS-C132931
Target:	WC1 Cell Surface Antigen
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Bovine
Immunogen:	Bovine lymphocytes.
Specificity:	Recognizes bovine WC1. Species cross-reactivity: sheep, goat.
Reactivity:	Bovine, Goat, Sheep
Purification:	Protein G purified
Presentation:	Supplied as a liquid in PBS, 1% BSA, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C. Sensitive to light.
Usage Summary:	Suitable for use in Flow cytometry: Neat-1:10; 10 ul labels 10 ⁶ cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences