

Human IgA Mouse anti-Human Monoclonal (RPE) (Hisa43) Antibody - LS-C134175 - LSBio	
CatalogID:	LS-C134175
Target:	Human IgA
Host	Human IgA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	Hisa43
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	Human IgA antibody was raised against Human
Immunogen:	Human IgA antibody was raised against human salivary proteins.
Specificity:	Cell surface expressed IgA.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	