

**CD14 Mouse anti-Human Monoclonal (FITC) (UCHM1) Antibody - LS-C134293 - LSBio**

<b>CatalogID:</b>	LS-C134293
<b>Target:</b>	CD14 molecule
<b>Synonyms:</b>	CD14 Antibody, CD14 molecule Antibody, CD14 antigen Antibody
<b>Host</b>	CD14 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	UCHM1
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD14 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD14 antibody was raised against human thymocytes/Sezary T cells.
<b>Specificity:</b>	Membrane bound CD14 molecule of about 53 kd. UCHM1 will also immunoprecipitate a soluble form of CD14.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	120 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences