


CD19 Mouse anti-Human Monoclonal (Biotin) (BU12) Antibody - LS-C134318 - LSBio

CatalogID:	LS-C134318
Target:	CD19 molecule
Synonyms:	CD19 Antibody, B-lymphocyte antigen CD19 Antibody, CD19 molecule Antibody, Differentiation antigen CD19 Antibody, T-cell surface antigen Leu-12 Antibody, B4 Antibody, CD19 antigen Antibody, CVID3 Antibody
Host	CD19 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	BU12
Conjugations:	Biotin
Immunogen Species:	CD19 antibody was raised against Human
Immunogen:	CD19 antibody was raised against human EB-4 Burkitt lymphoma cell line.
Specificity:	CD19 molecule of about 95 kd.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. DO NOT FREEZE.
Usage Summary:	Function: Induces adhesion of B cells.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg

Flow Cytometry Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences