

CD22 Mouse anti-Human Monoclonal (FITC) (RFB4) Antibody - LS-C134334 - LSBio

CatalogID:	LS-C134334
Target:	CD22 molecule
Synonyms:	CD22 Antibody, CD22 molecule Antibody, CD22 antigen Antibody, SIGLEC2 Antibody, SIGLEC-2 Antibody, T-cell surface antigen Leu-14 Antibody, B-cell receptor CD22 Antibody, BL-CAM Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD22 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	RFB4
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD22 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD22 antibody was raised against human tonsil lymphocytes.
Specificity:	CD22 molecule of about 140 kd.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	IHC - Paraffin, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences