

FCER2 / CD23 Mouse anti-Human Monoclonal (RPE) (Bu38) Antibody - LS-C134340 - LSBio	
CatalogID:	LS-C134340
Target:	Fc fragment of IgE, low affinity II, receptor for (CD23) (FCER2)
Synonyms:	FCER2 Antibody, CD23 antigen Antibody, CD23A Antibody, CLEC4J Antibody, CD23 Antibody, Fc-epsilon-RII Antibody, IGEBF Antibody, Lymphocyte IgE receptor Antibody, BLAST-2 Antibody, FCE2 Antibody
Host	FCER2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Bu38
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	FCER2 / CD23 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	FCER2 / CD23 antibody was raised against human peripheral myeloma cells.
Specificity:	CD23 molecule of 45 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	IHC - Paraffin, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	