


PECAM-1 / CD31 Mouse anti-Human Monoclonal (FITC) (158-2B3) Antibody - LS-C134380 - LSBio	
CatalogID:	LS-C134380
Target:	platelet/endothelial cell adhesion molecule 1 (PECAM1)
Synonyms:	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
Host	PECAM1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	158-2B3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PECAM-1 / CD31 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	PECAM-1 / CD31 antibody was raised against stimulated human leukocytes.
Specificity:	Domain 1 of CD31.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Blocks homophilic interaction and heterophilic transendothelial migration.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Flow cytometry of PECAM / CD31 antibody

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences