

PECAM-1 / CD31 Mouse anti-Human Monoclonal (RPE) (158-2B3) Antibody - LS-C134381 - LSBio	
<b>CatalogID:</b>	LS-C134381
<b>Target:</b>	platelet/endothelial cell adhesion molecule 1 (PECAM1)
<b>Synonyms:</b>	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
<b>Host</b>	PECAM1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	158-2B3
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	PECAM-1 / CD31 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	PECAM-1 / CD31 antibody was raised against stimulated human leukocytes.
<b>Specificity:</b>	Domain 1 of CD31.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Gel filtration
<b>Presentation:</b>	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
<b>Usage Summary:</b>	Function: Blocks homophilic interaction and heterophilic transendothelial migration.
<b>Uses:</b>	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	120 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	