


CD34 Mouse anti-Human Monoclonal (APC) (43A1) Antibody - LS-C134398 - LSBio

CatalogID:	LS-C134398
Target:	CD34 molecule
Synonyms:	CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody
Host	CD34 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG3,k
Clone Name:	43A1
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	CD34 antibody was raised against Human
Immunogen:	CD34 antibody was raised against human KG-1a tumor cells.
Specificity:	Class II epitope of the 110 kd CD34 molecule which is resistant to neuraminidase and sensitive to glycoprotease.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	IHC - Paraffin, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Flow cytometry of CD34 antibody

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences