


CD37 Mouse anti-Human Monoclonal (RPE) (IPO-24) Antibody - LS-C134412 - LSBio

CatalogID:	LS-C134412
Target:	CD37 molecule
Synonyms:	CD37 Antibody, CD37 antigen Antibody, CD37 molecule Antibody, Leukocyte surface antigen CD37 Antibody, gp52-40 Antibody, TSPAN26 Antibody, Leukocyte antigen CD37 Antibody, Tetraspanin-26 Antibody, Tspan-26 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD37 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	IPO-24
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD37 antibody was raised against Human
Immunogen:	CD37 antibody was raised against human spleen cells from a patient with hairy cell leukemia.
Specificity:	IPO-24 recognizes the CD37 molecule of about 42-52 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Flow cytometry of CD37 antibody

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences