

ITGA2B / CD41 Mouse anti-Human Monoclonal (RPE) (96.2C1) Antibody - LS-C134432 - LSBio	
CatalogID:	LS-C134432
Target:	integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41) (ITGA2B)
Synonyms:	ITGA2B Antibody, BDPLT2 Antibody, CD41B Antibody, GT Antibody, gp2B Antibody, Integrin alpha-IIb Antibody, ITGAB Antibody, GPIIb Antibody, HPA3 Antibody, Platelet-specific antigen BAK Antibody, CD41 Antibody, CD41 antigen Antibody, GPalpha IIb Antibody, GTA Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGA2B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	96.2C1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	ITGA2B / CD41 antibody was raised against Human
Immunogen:	ITGA2B / CD41 antibody was raised against human PBL.
Specificity:	CD41a of about 140 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/25/2014 © 2014 LifeSpan BioSciences	