

CD43 Mouse anti-Human Monoclonal (RPE) (DFT-1) Antibody - LS-C134437 - LSBio

CatalogID:	LS-C134437
Target:	sialophorin (SPN)
Synonyms:	SPN Antibody, GPL115 Antibody, Leukosialin Antibody, Leukocyte sialoglycoprotein Antibody, Galactoglycoprotein Antibody, GALGP Antibody, Sialophorin Antibody, LSN Antibody, CD43 Antibody, CD43 antigen Antibody
Host	SPN antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	DFT-1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD43 antibody was raised against Human
Immunogen:	CD43 antibody was raised against human KG-1 tumor cells.
Specificity:	Epitope A of the CD43 molecule on pMN of about 135 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Partially induces apoptosis in hemopoietic progenitor cells and also induces hemopoietic aggregation.
Uses:	IHC, Western blot, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences