

**CD45R Mouse anti-Human Monoclonal (Azide-free) (351C5) Antibody - LS-C134452 - LSBio**

<b>CatalogID:</b>	LS-C134452
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgM
<b>Clone Name:</b>	351C5
<b>Modifications:</b>	Azide-free
<b>Immunogen Species:</b>	CD45R antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45R antibody was raised against peripheral blood mononuclear cells.
<b>Specificity:</b>	CD45R isoforms of about 200-220 kd.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Gel filtration
<b>Presentation:</b>	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Open under aseptic conditions.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences