

CD45 Mouse anti-Human Monoclonal (FITC) (C11) Antibody - LS-C134455 - LSBio

CatalogID:	LS-C134455
Target:	protein tyrosine phosphatase, receptor type, C (PTPRC)
Synonyms:	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
Family / Subfamily:	Protein Phosphatase / R1/R6
Host	PTPRC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	C11
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45 antibody was raised against peripheral blood mononuclear cells.
Specificity:	All isoforms of CD45.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences